stratone

User's Manual
ENGLISH

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M-70-15-0862(02)

Promedon

USER MANUAL

DESCRIPTION

¡PLEASE NOTE! This product is packed sterile (ETO sterilization) and intended for disposal after use. The product must be stored in a cool dry place. It should not be used if it is damaged. Only a doctor may use this product.

DESCRIPTION AND STRUCTURE

The universal stone basket is made of nitinol and consists of three main components:

1. Stone basket.
2. Sheath.
3. Detachable handle with screw.

INDICATIONS

Endoscopic removal of urological stones.

POSSIBLE COMPLICATIONS

- Oedema.
- Dissection.
- Urosepsis.
- Perforation of the ureter.
- Latrogenic related infections.

PREPARATION

- Check packaging for damage.
- Remove product from packaging.
- Check product for damage.

HANDLING

- Introduce the optical instrument into the urinary tract.
- By means of the optical instrument, and with closed basket the Universal Stone Basket is manoeuvred so that the tip of the basket lies just above the stone.
- By means of pushing the slide of the handle the sheath is pulled back, thus opening the basket.
- Manipulate the basket until the stone is caught inside the basket.
- As soon as the stone is in the basket the sheath is moved forward by means of pulling the slide of the handle (closing the basket) until the stone is fixed.
- Small stones which do not damage the walls of the urinary tract can now be removed.
- Larger stones which would damage the walls of the urinary tract should be removed by some other method.

NOTE: after detachment of the handle, the universal stone basket furthermore permits easy removal of the sheath through the optical instrument, whereby the basket may simply remain in place (See Detachment of the handle).

Detachment of the handle: Unscrew the screw in the handle and remove the sheath releasing the Luer connector (red) attached to the handle.

SYMBOLS USED IN THE LABELS

CATALOGUE NUMBER

BATCH NUMBER

SEE INSTRUCTIONS FOR USE

EXPIRATION DATE

DO NOT REUSE

DATE OF MANUFACTURE

STERILE: STERILIZATION METHOD
ETHYLENE OXIDE

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